ARGENTA

## *Standard Q*

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### CRP


STANDARD™ F CRP

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

REF<sub>F-CRP</sub>

STANDARD™

Zakres:  
WB:1-150 mg/L  
S/P : 1-130

STANDARD F

### PCT FIA

STANDARD™ F PCT FIA

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

REF<sub>F-PCT-02</sub>

STANDARD™

WB/S/P: 0,05-50 ng/ml

STANDARD F

### PCT FIA (Serum)


STANDARD™ F PCT FIA (Serum)

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

REF<sub>F-PCT-01</sub>

STANDARD™

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Źródło: Protein Data Bank

ARGENTA

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### STANDARD F β-hCG FIA

STANDARD™ F β-hCG FIA

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST


REF F-HCG

STANDARD™

### STANDARD F LH FIA

STANDARD™ F LH FIA

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

Źródło: Protein Data Bank

### STANDARD F TSH FIA

STANDARD™ F TSH FIA

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

REF F-TSH

STANDARD™

### STANDARD F ft4

STANDARD™ F ft4

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

REF F-ft4-01

STANDARD™

### STANDARD F T4

STANDARD™ F T4

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE YOU PERFORM THE TEST

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
Group	Viral Transport Medium
VTM – Group 1	COPAN UTM™
VTM – Group 2	STANDARD™ Transport Medium
	BD Universal Viral Transport

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
## STANDARD F COVID-19 Ag FIA | Test Procedure


## ✓ Test Procedure – STANDARD Test Mode

1 Collect specimen

2 Insert test device

3 Mix specimen into extraction buffer tube

4 Select specimen type


5 Apply specimen and push START


6 Test result will be displayed in 30 minutes


1) Excessive amount of specimen

2) Adequate amount of specimen

3) Specimen collection technique when patient is undergoing excessive nasal discharge in their nasal cavity.

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# Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays

Interim guidance

11 September 2020

1. SARS-CoV-2 Ag-RDTs that meet the minimum performance requirements of  $\geq 80\%$  sensitivity and  $\geq 97\%$  specificity compared to a NAAT reference assay<sup>1</sup> can be used to diagnose SARS-CoV-2 infection in a range of settings where NAAT is unavailable or where prolonged turnaround times preclude clinical utility.

To optimize performance, testing with Ag-RDTs should be conducted by trained operators in strict accordance with the manufacturer's instructions and within the first 5-7 days following the onset of symptoms.

Although more evidence is needed on real-world performance and operational aspects, Ag-RDTs are most likely to perform well in patients with high viral loads (Ct values  $\leq 25$  or  $> 10^6$  genomic virus copies/mL) which usually appear in the pre-symptomatic (1-3 days before symptom onset) and early symptomatic phases of the illness (within the first 5-7 days of illness) (14, 15). This

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Specimen	Concentration	Result analysis
Positive Nasopharyngeal swab specimen	1 x Limit of Detection: 7.81 X 10 <sup>1.2</sup> TCID <sub>50</sub> /ml	100% Sensitivity (30/30)
	2x Limit of Detection: 1.56 X 10 <sup>2.2</sup> TCID <sub>50</sub> /ml	
	4x Limit of Detection: 3.13 X 10 <sup>3.2</sup> TCID <sub>50</sub> /ml	
Negative Nasopharyngeal swab specimen	N/A	100% Specificity (30/30)

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## Clinical Study : Korea

- Product : STANDARD F COVID-19 Ag FIA
- A total of 155(55 Positive, 100 Negative) specimens were tested
- Sensitivity: **89.09%**
- Specificity: **96%**
- Test Site: YUMC

		Real-time PCR		Total
		Positive	Negative	
STANDARD F COVID-19 Ag Test	Positive	49	4	53
	Negative	6	96	102
Total		55	100	155

### Conclusion


In conclusion, through this clinical performance evaluation, it was confirmed that “STANDARD F COVID-19 Ag FIA”, has high clinical sensitivity and specificity compared to the evaluation criteria, and the correlation with the symptom onset date was confirmed. When confirmed, it was confirmed that the faster the test was performed at the onset of symptoms, the higher the test performance was. Therefore, the “STANDARD Q COVID-19 Ag Test”, a medical device for clinical performance testing, will be very useful as a diagnostic product that quickly and screens for novel coronavirus infection-19 (COVID-19).

## Clinical Study : Brazil

### < STANDARD F COVID-19 Ag FIA >

- A total of 24 (5 Positive, 19 Negative) specimens were tested
- Sensitivity: **91.70% (CI = 95%)**
- Specificity: **100.00% (CI = 95%)**
- Test Site: **PPC/UERJ**

NP swab in VTM		STANDARD F COVID-19 Ag FIA		Total
		Positive	Negative	
SARS-CoV-2 rt-PCR	Positive	4	1	5
	Negative	0	19	19
<b>Total</b>		<b>4</b>	<b>20</b>	<b>24</b>
<b>Clinical Sensitivity</b>		<b>91.70% (CI = 95%)</b>		
<b>Clinical Specificity</b>		<b>100% (CI = 95%)</b>		

#### Conclusion

As shown above, we can conclude that the Covid-19 Ag FIA Test kit, manufactured by SD Biosensor, is suitable for performing the rapid antigen detection test for Covid-19.

## Clinical Study :Spain

### < STANDARD F COVID-19 Ag FIA >

- A total of 665 (111 Positive - Ct≤30, 14 Positive Ct>30, 538 Negative) specimens were tested
- Sensitivity (Ct≤30) : **93.70% (95% CI: 87.44 – 97.43%)**
- Specificity: **99.63 % (95% CI: 98.50% - 99.96%)**
- Test Site: **University Hospital Marques de Valdecilla**

NP swab in UTM		STANDARD F COVID-19 Ag FIA		Total
		Positive	Negative	
SARS-CoV-2 rt-PCR	Positive	104	7	111
	Negative	2	536	538
Total		106	543	649
Clinical Sensitivity		93.70% (95% CI: 87.44 – 97.43%)		
Clinical Specificity		99.63% (95% CI: 98.50% - 99.96%)		

RT-PCR Ct value	No of patients	Sensitivity(%)
Ct < 15	6	100% (6/6)
Ct 16 – 20	41	100% (41/41)
Ct 21 – 25	47	93.62% (44/47)
Ct 26 – 30	17	76.5% (13/17)
Clinical Sensitivity (Ct≤30)		<b>93.70% (104/111)</b>

- The positive predictive value (PPV) was 98.1% (92.7, 99.7).
- The negative predictive value (NPV) was 96.4% (94.4, 97.7).
- The sensitivity including 14 patients with RT-PCR Ct value >30 was 84.00%.

## Clinical Study :India

### < STANDARD F COVID-19 Ag FIA >


- A total of 303 (87 Positive, 216 Negative) specimens were tested
- Sensitivity: **93.10% (95% CI: 85.59 – 97.43%)**
- Specificity: **99.07 % (95% CI: 96.70% - 99.89%)**
- Test Site: **All India Institute of Medical Sciences**

NP swab in UTM		STANDARD F COVID-19 Ag FIA		Total
		Positive	Negative	
SARS-CoV-2 rt-PCR	Positive	81	6	87
	Negative	2	214	216
Total		83	220	303
Clinical Sensitivity		93.10% (95% CI: 85.59 – 97.43%)		
Clinical Specificity		99.07% (95% CI: 96.70% - 99.89%)		

*Clinical Study : U.A.E. DUBAI*< STANDARD F COVID-19 Ag FIA >

- A total of 63( 32 Positive, 31 Negative) specimens were tested
- Sensitivity: **93.75%** (95% CI: 80.32 – 99.28%)
- Specificity: **96.77%** (95% CI: 83.30% - 99.92%)
- Test Site: **RAK hospital Laboratory**

NP swab in VTM		STANDARD F COVID-19 Ag FIA		Total
		Positive	Negative	
SARS-CoV-2 RT-PCR	Positive	30	2	32
	Negative	1	30	31
Total		31	32	63

# Summary of WHO Interim Guidance (11 Sept, 2020)

[Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays]

2020.12  
F-COVID-WHO\_R.02

## STANDARD F COVID-19 Ag FIA

### 1. General Recommendation


No.	WHO recommendations	STANDARD F COVID-19 Ag FIA
1	SARS-CoV-2 Ag RDTs that meet the minimum performance requirements of $\geq 80\%$ sensitivity and $\geq 97\%$ specificity compared to a NAAT reference assay.	According to clinical evaluations conducted in nine different countries around the world, total clinical sensitivity & specificity of STANDARD F COVID-19 Ag FIA is <b>94.09%</b> and <b>98.52%</b> respectively.

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





Result	COI (Cutoff index) value	SARS-COV-2 Ag
Positive	COI $\geq$ 1.0	Positive for SARS-COV-2 Ag
Negative	COI < 1.0	Negative for SARS-COV-2 Ag
Invalid	COI value is not displayed	Retest should be performed with a new test device and a new patient's specimen.

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Mobile

Personal computer

STANDARD  
The Global Standard Founder

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