

A microscopic view of foodborne pathogens, showing several rod-shaped bacteria with long, thin flagella. The bacteria are rendered in a light blue/white color against a dark background. A horizontal orange band is overlaid across the middle of the image, containing the title text.

FOODBORNE PATHOGENS


Eurofins Technologies - @Argenta

05.10.2022

Ricardo Silva e Sousa
RicardoSilvaeSousa@eurofins.de
Jr. Product Manager

- Why are we here? Statistics on Salmonella, Listeria & STEC E. coli
- Foodborne pathogens – regulation, methods & markets
- Portfolio overview – What is the *BACGene* product line
- Extra troubleshooting methods
- Available tools
- Technical overview

Why are we here?




Salmonella



Technologies

	2019	2018	2017	2016	2015	Data source
Humans						
Total number of confirmed cases	87,923	91,858	91,587	94,425	94,477	ECDC
Total number of confirmed cases/100,000 population (notification rates)	20.0	20.1	19.7	20.5	21.0	ECDC
Number of reporting MS	28	28	28	28	28	ECDC
Infection acquired in the EU	58,271	59,763	59,642	52,852	51,898	ECDC
Infection acquired outside the EU	6,343	6,376	6,001	6,466	6,830	ECDC
Unknown travel status or unknown country of infection	23,309	25,719	25,944	35,107	35,749	ECDC
Number of outbreak-related cases	9,169	11,631	9,607	11,428	8,531	EFSA
Total number of outbreaks	926	1,588	1,241	1,372	1,216	EFSA
Food						
Meat and meat products						
Number of sampling units	525,704	433,197	380,000	285,564	211,072	EFSA
Number of reporting countries	28	28	28	27	27	EFSA
Milk and milk products						
Number of sampling units	46,797	44,078	30,796	24,337	29,034	EFSA
Number of reporting countries	25	24	24	24	22	EFSA
Fish and fishery products						
Number of sampling units	14,010	17,123	13,507	12,287	11,373	EFSA
Number of reporting countries	24	22	22	21	22	EFSA
Eggs and egg products						
Number of sampling units	12,093	10,611	15,435	10,933	9,650	EFSA
Number of reporting countries	21	21	23	20	19	EFSA
Fruits and vegetables (and juices)						
Number of sampling units	17,068	10,888	7,579	7,515	6,797	EFSA
Number of reporting countries	22	22	25	20	22	EFSA



Summary of Salmonella monitoring results, by major meat and meat products categories and by sampling stage in the food chain, EU, 2019

Summary of Salmonella statistics related to humans and major food categories, EU, 2015–2019

Salmonella:

- Second most common zoonosis in humans in the EU
- poultry meats – where *Salmonella* was most frequently reported
- egg and eggs products ranked first of food vehicles
- Previous decreasing trend has stabilized
- The fight is not being won. Rates of cases are flat.

Country	National coverage ^(a)	Data format ^(a)	2019		2018		2017		2016		2015	
			Total cases	Confirmed cases & rates	Confirmed cases & rates	Confirmed cases & rates	Confirmed cases & rates	Confirmed cases & rates				
									Cases	Rate	Cases	Rate
Poland	Y	C	8,919	8,373 22.0	9,064 23.9	8,921 23.5	9,718 25.6	8,245 21.7				

Reported human cases of salmonellosis and notification rates per 100,000 population in the EU/EFTA, by country and year, 2015–2019

EFSA Journal 2021;19(2):6406

Why are we here?



Listeria



Technologies

	2019	2018	2017	2016	2015	Data source
Humans						
Total number of confirmed cases	2,621	2,545	2,475	2,500	2,183	ECDC
Total number of confirmed cases/100,000 population (notification rates)	0.46	0.47	0.48	0.47	0.43	ECDC
Number of reporting MS	28	28	28	28	28	ECDC
Infection acquired in the EU	1,817	1,640	1,639	1,539	1,450	ECDC
Infection acquired outside the EU	12	8	4	6	7	ECDC
Unknown travel status or unknown country of infection	792	897	832	955	726	ECDC
Number of outbreak-related cases	349	159	39	27	233	ECDC
Total number of outbreaks	21	14	10	6	15	EFSA
RTE food categories^(a)						
RTE milk and milk products	N = 62,019; 23 MS	N = 59,313; 23 MS	N = 56,428; 25 MS	N = 34,850; 26 MS	N = 45,996; 24 MS	EFSA
RTE meat and meat products	N = 64,666; 22 MS	N = 57,861; 22 MS	N = 45,219; 24 MS	N = 25,195; 21 MS	N = 25,396; 22 MS	EFSA
RTE fish and fishery products	N = 13,376; 22 MS	N = 14,081; 22 MS	N = 12,604; 24 MS	N = 6,601; 21 MS	N = 7,986; 25 MS	EFSA
Other RTE food products	N = 76,657; 24 MS	N = 25,179; 22 MS	N = 23,915; 23 MS	N = 21,085; 22 MS	N = 25,544; 23 MS	EFSA
RTE foods intended for infants and for special medical purposes	N = 1,721; 18 MS	N = 1,663; 18 MS	N = 1,462; 20 MS	N = 1,274; 16 MS	N = 1,754; 12 MS	EFSA

Summary statistics on human invasive *L. monocytogenes* infections and on sampled major RTE food categories in the EU, 2015–2019

Country	2019						2018		2017		2016		2015	
	National coverage ^(a)	Data format ^(a)	Total cases	Confirmed cases and rates		Confirmed cases and rates		Confirmed cases and rates		Confirmed cases and rates		Confirmed cases and rates		
				Cases	Rate	Cases	Rate	Cases	Rate	Cases	Rate	Cases	Rate	
Poland	Y	C	121	121	0.32	128	0.34	116	0.31	101	0.27	70	0.18	

Reported cases of human invasive listeriosis and notification rates per 100,000 population in the EU/EFTA, by country and year, 2015–2019

Listeria:

- The EU trend of confirmed listeriosis cases remained stable (flat) in 2015–2019 after a long period of an increasing trend.
- The overall EU case fatality was high (17.6%) - listeriosis one of the most serious food-borne diseases under EU surveillance

The battlefield is currently on the Ready To Eat products.

There was an increase of:

- 12% of the sampling units tested for 'RTE meat and meat products'
- 204% for 'other RTE food products'.
 - higher number of samples were tested for 'bakery products' (+75%),
 - 'broiler meat and meat products thereof' (+304%)
 - fruit and vegetables (+79%).

Why are we here?



STEC *E. coli*



Technologies

Reported human cases of STEC infections and notification rates per 100,000 population in the EU/EFTA, by country and year, 2015–2019

	2019	2018	2017	2016	2015	Data source
Humans						
Total number of confirmed cases	7,775	8,161	5,958	6,474	5,929	ECDC
Total number of confirmed cases/100,000 population (notification rates)	2.21	2.28	1.67	1.79	1.65	ECDC
Number of reporting MS	27	28	28	28	28	ECDC
Infection acquired in the EU	4,835	5,783	4,747	4,037	3,991	ECDC
Infection acquired outside the EU	750	693	525	339	532	ECDC
Unknown travel status or unknown country of infection	2,190	1,685	686	2,098	1,406	ECDC
Number of food-borne outbreak-related cases	273	390	260	737	676	EFSA
Total number of food-borne outbreaks	42	50	48	43	70	EFSA
Food						
All						
Number of sampling units	25,030	20,498	19,351	17,977	13,777	EFSA
Number of reporting MS	22	20	22	17	17	EFSA
Meat and meat products						
Number of sampling units	14,110	9,250	10,706	8,771	7,865	EFSA
Number of reporting MS	20	17	18	17	15	EFSA
Milk and milk products						
Number of sampling units	5,479	5,339	3,485	3,773	3,005	EFSA
Number of reporting MS	13	14	10	11	8	EFSA
Fruits and vegetables (and juices)						
Number of sampling units	2,658	3,339	2,295	1,475	1,384	EFSA
Number of reporting MS	13	13	15	11	10	EFSA

Summary of STEC statistics related to humans and to major food categories, EU, 2015–2019

Country	2019			2018		2017		2016		2015			
	National coverage ^(a)	Data format ^(a)	Total cases	Confirmed cases and rates		Confirmed cases and rates		Confirmed cases and rates		Confirmed cases and rates			
				Cases	Rate	Cases	Rate	Cases	Rate	Cases	Rate		
Poland	Y	C	17	14	0.04	6	0.01	4	0.01	4	0.01	0	0.00

STEC *E. coli*:

- The EU/EEA trend has been increasing from 2015 to 2019.
- STEC - 3rd most frequent bacterial agent detected in food-borne outbreaks in the EU:
 - 42 outbreaks, 273 cases, 50 hospitalisations and 1 death (2019)
- Sprouted - tested by six MS with no positive STEC results from 331 official samples.
 - An EU regulation has been in force since 2013.
- Overall, STEC was most commonly found in:
 - meat of different types derived from different animal species (4.1% STEC-positive),
 - followed by 'milk and dairy products' (2.1%)
 - while 'fruits and vegetables' was the least contaminated category (0.1%).


Food category	Micro-organisms/their toxins, metabolites
1.1. Ready-to-eat foods intended for infants and ready-to-eat foods for special medical purposes (*)	<i>Listeria monocytogenes</i>
1.2. Ready-to-eat foods able to support the growth of <i>L. monocytogenes</i> , other than those intended for infants and for special medical purposes	<i>Listeria monocytogenes</i>
1.3. Ready-to-eat foods unable to support the growth of <i>L. monocytogenes</i> , other than those intended for infants and for special medical purposes (*) (†)	<i>Listeria monocytogenes</i>
1.4. Minced meat and meat preparations intended to be eaten raw	<i>Salmonella</i>
1.5. Minced meat and meat preparations made from poultry meat intended to be eaten cooked	<i>Salmonella</i>
1.6. Minced meat and meat preparations made from other species than poultry intended to be eaten cooked	<i>Salmonella</i>
1.7. Mechanically separated meat (MSM) (‡)	<i>Salmonella</i>
1.8. Meat products intended to be eaten raw, excluding products where the manufacturing process or the composition of the product will eliminate the salmonella risk	<i>Salmonella</i>

Food category	Micro-organisms/their toxins, metabolites
1.9. Meat products made from poultry meat intended to be eaten cooked	<i>Salmonella</i>
1.10. Gelatine and collagen	<i>Salmonella</i>
1.11. Cheeses, butter and cream made from raw milk or milk that has undergone a lower heat treatment than pasteurisation (10)	<i>Salmonella</i>
1.12. Milk powder and whey powder (10)	<i>Salmonella</i>
1.13. Ice cream (11), excluding products where the manufacturing process or the composition of the product will eliminate the salmonella risk	<i>Salmonella</i>
1.14. Egg products, excluding products where the manufacturing process or the composition of the product will eliminate the salmonella risk	<i>Salmonella</i>
1.15. Ready-to-eat foods containing raw egg, excluding products where the manufacturing process or the composition of the product will eliminate the salmonella risk	<i>Salmonella</i>
1.16. Cooked crustaceans and molluscan shellfish	<i>Salmonella</i>
1.17. Live bivalve molluscs and live echinoderms, tunicates and gastropods	<i>Salmonella</i>


Food category	Micro-organisms/their toxins, metabolites
1.18. Sprouted seeds (ready-to-eat) (12)	<i>Salmonella</i>
1.19. Pre-cut fruit and vegetables (ready-to-eat)	<i>Salmonella</i>
1.20. Unpasteurised fruit and vegetable juices (ready-to-eat)	<i>Salmonella</i>
1.21. Cheeses, milk powder and whey powder as referred to in the coagulase-positive staphylococci criteria in Chapter 2.2 of this Annex	Staphylococcal enterotoxins
1.22. Dried infant formulae and dried dietary foods for special medical purposes intended for infants below six months of age, as referred to in the Enterobacteriaceae criterion in Chapter 2.2 of this Annex	<i>Salmonella</i>
1.23. Dried infant formulae and dried dietary foods for special medical purposes intended for infants below six months of age, as referred to in the Enterobacteriaceae criterion in Chapter 2.2 of this Annex	<i>Enterobacter sakazakii</i>
1.24. Live bivalve molluscs and live echinoderms, tunicates and gastropods	<i>E.coli</i> (14)
1.25. Fishery products from fish species associated with a high amount of histidine (15)	Histamine

Regulation:


- EC no 2073/2005 is the most relevant document towards EU regulation
- Most of the regulation requires absence in 10/25g of either *Salmonella* or *Listeria monocytogenes*.
- There is almost no EU regulation for STEC (addendum for sprouts since 2013)
- *Campylobacter* also added for broiler carcasses (Chapter 2 of Annex I)
- Other testing results from local regulations or market pressure – **bottom up approach**.



Salmonella




Listeria




*STEC
E coli*

Background Information: *Salmonella*

- Gram negative
- *Enterobacteriaceae*
- 2 species: *S. enterica* & *S. bongori*
- Six *S. enterica* sub-species



Salmonella




Salmonella

- Symptoms include diarrhea, fever and stomach cramps
- Symptoms start within 6 - 72 hours and last between 4 - 7 days

Background Information: *Listeria*

- Gram positive
- Belonging to the phylum Firmicutes
- Resistant to stress (temp, salt, pH)



Known pathogenic species:
L. monocytogenes (13 serovars)

Listeria monocytogenes

- Listeriosis is a serious infection (fever – flu-like symptoms)
- Most likely to sicken pregnant women and their newborns, elderly and immunocompromised people

- Gram negative
- *Enterobacteriaceae*
- Common gut flora
- Most non-pathogenic




O-Group

O157
O145
O111
O121
O103
O45
O26

Top 7


O157:H7



Shiga toxins
(*stx1* & *stx2*)




Intimin
(*eae*)




STEC – Shiga toxin-producing *E. coli*

- Symptoms: bloody diarrhoea, cramps, vomiting
- Severe diseases such as haemolytic uremic syndrome (HUS) or haemorrhagic colitis (HC)
- Highly infectious (infection dose: <10 cells)




Target pathogens:

- Majority of testing target towards *Salmonella* and *Listeria*.
- Despite *Campylobacter* being the most common cause for gastroenteritis in the EU is one of the least tested.
- Honorable mention to *Cronobacter* which is very relevant for the infant formula industry



- Protein segment includes
 - Meat, fish, poultry and eggs
- Processed food highest variety
- Fruits/vegetables ranges third
- Lowest dairy



Testing methods:

- Molecular methods with high growth potential by entering other markets.
- Company size/volumes must reach a threshold to consider an in-house laboratory.
- Chromogenic is sizeable in EU. Highly used in the UK market which is very price sensitive.

Immunoassay = ELISA-based Antibody systems


Molecular = real-time PCR

Traditional media = plate counts

Chromogenic = targeted visualization of microcolonies

Traditional media vs ELISA vs qPCR

	ELISA	ISO	PCR	
	Immunoassay	Culture	Molecular	
Ease-of-use	++	-	+	Ease-of-use
Technical skills	-	-*	+	Technical skills
Costs	\$	\$	\$\$	Costs
Time to Neg. Result	+	-	++	Time to Neg. Result
Specificity	+/-	+*	++	Specificity



Testing methods:

In the end not every customer can use molecular methods.

There is a market and a place for every method.

- Traditional media
 - Inexpensive
 - Long TTR
 - Requires a trained eye from an experienced technician.
- ELISA
 - Easy to use and interpret especially with the automated methods available (robots and readers)
 - Suboptimal specificity
- qPCR
 - Fastest by a large margin
 - Increasingly simple
 - Very high specificity

PCR **BACGene** *Overview*



BACGene: Real-time PCR based detection of pathogens in food, feed and environmental samples

Short Time to Results

- **Short** enrichment & **Absence** of secondary enrichment step for standard matrices
- **Simple** DNA preparation

High Reliability

- Probe based assay => excellent **specificity**
- Positive and negative PCR controls => PCR target detection and **absence of contamination**
- **Internal Positive Control (IPC)** for verification of PCR performance
- **PREraser** – certified free-DNA removal & **UNG** – cross contamination prevention (STEC only)

Proven performance

- Certified by **independent** organisations (AFNOR, AOAC)
- **Automated** data analysis (Excel Evaluation Sheet)

High Flexibility

- **Open system**: certified on AriaMX and Bio-Rad CFX96
- One PCR program for **all** BACGene kits


Overview BACGene portfolio

BACGene			Extras
Salmonella spp. <i>Salmonella spp.</i>	Listeria spp. <i>Listeria spp.</i>	E. coli O157:H7 <i>E. coli</i> serotype O157:H7	PREraser BACGene (Elimination Free-DNA)
	Listeria monocytogenes <i>Listeria monocytogenes</i>	E. coli STEC Multiplex STEC virulence genes stx1, stx2 & eae	
	Listeria Multiplex <i>Listeria monocytogenes</i> and <i>Listeria spp.</i>	E. coli Serotype 1 STEC serotypes O111, O103 & O26	Salmonella Tranoroa (Cross Contamination Troubleshooting)
	E. coli Serotype 2 STEC serotypes O145, O121 & O45		




The BACGene Solution

Enrichment



Lysis




PCR



Thermocycler



Readout



Sample		Reference		Result	
ID	Barcode	ID	Barcode	Value	Unit
101	101	101	101	20.00	2.000 mg/dl
102	102	102	102	20.00	2.000 mg/dl
103	103	103	103	20.00	2.000 mg/dl
104	104	104	104	20.00	2.000 mg/dl
105	105	105	105	20.00	2.000 mg/dl
106	106	106	106	20.00	2.000 mg/dl
107	107	107	107	20.00	2.000 mg/dl
108	108	108	108	20.00	2.000 mg/dl
109	109	109	109	20.00	2.000 mg/dl
110	110	110	110	20.00	2.000 mg/dl

BACGene kits

Evaluation sheets

Media Development
Ongoing

Agreement
Agilent (Aria MX)

	Salmonella	Listeria	E.coli
Media	BPW	Actero	mTSB
Sample size	25 g 375 g	25 g	25 g 375 g
Dilutions	1:10	1:10	1:3 (meat) 1:5 (leafy greens)
Subbing	Chocolate (NF milk, not BPW) -> BHI		

Enrichment – What is our Matrix Scope?

BACGene *Salmonella* spp. - AFNOR Certified Protocols

Sample Matrix	Milk and dairy (25 g)	All foods (25 g)	Feed (25 g)	Pet food (up to 375 g)	Milk powders & infant formula (up to 375 g)	Milk powders & infant formula + probiotics (up to 375 g)	Environmentals
	Media	1:10 BPW			1:10 BPW pre-warmed		1:10 BPW ² pre-warmed
Enrichment T °C	41,5 ± 1°C		37 ± 1°C				
Time	18 - 24 h	16 - 24 h	18 - 24 h				



BACGene *Listeria* spp., *monocytogenes* & Multiplex - AFNOR Certified Protocols

Sample Matrix	All foods (25 g)	Dust & process water (25 g/mL)	Environmentals
	Media	1:10 Actero™ pre-warmed	
Enrichment T °C	37 ± 1°C		
Time	18 - 24 h		



	Salmonella	Listeria	E.coli
Lysis Buffer	Lysis buffer I	Lysis buffer L	Lysis buffer S
Enzyme	Proteinase K	Proteinase K + Lysozyme	Proteinase K
Incubation	37±2°C - 20 min + 95±5°C - 10 min		
Equipment needed	Heating block		

PCR & Thermocyclers

	Salmonella	Listeria	E.coli
Run time	~ 1:30 h:min		
Thermal profile	Same		
Thermo Cycler	Agilent AriaMX™ Bio-rad CFX 96 touch™ Standard Bio-rad CFX 96 touch™ Deep Well		
Certification	Afnor & AOAC	Afnor & AOAC	Afnor (O157:H7) & AOAC



Data Interpretation

Example - BACGene Salmonella spp.

<i>Salmonella</i>	Internal Positive Control	Final results	
Reaction positive	Valid/Invalid (<i>Not significant</i>)	Positive	123 +
Reaction negative	Valid	Negative	456 -
Reaction negative	Invalid	Questionable*	789 Q

* Refer to troubleshooting section in the *BACGene Salmonella spp.* manual

PCR **BACGene**
troubleshooting

Questions:

- Does a positive result comes from a cross contamination?
 - **BACGene Salmonella Tranoroa**
- Does a positive result comes from viable or non-viable pathogen?
 - **PREraser BACGene**

Questions:

- Does a positive result comes from a cross contamination?
 - **BACGene Salmonella Tranoroa**
- Does a positive result comes from viable or non-viable pathogen?
 - *PREraser BACGene*



Ideal routine positive control strain

- Not naturally occurring in food matrices
- Easily and quickly discriminated from other serotypes
- Same growth conditions as food related serotypes



***Salmonella enterica* serovar Tranoroa**


BACGene *Salmonella* Tranoroa a rapid test from the enrichment broth

Can be run after BACGene *Salmonella* spp. from the same lysate to **confirm absence of S. Tranoroa**

Evaluation sheets BACGene includes *Salmonella* Tranoroa



Integration of BACGene *Salmonella* Tranoroa (enrichment positive control)




Questions:

- Does a positive result comes from a cross contamination?
 - *BACGene Salmonella Tranoroa*
- Does a positive result comes from viable or non-viable pathogen?
 - ***PREraser BACGene***


DNA from Non Viable Microorganisms

PCR in a Microbiological Lab:


Microorganism



Microorganism & DNA




DNA





- Enzymatic pre-treatment prior to lysis
- In case of high rate of not confirmable PCR positive results due to free DNA



PCR **BAC** *Gene*
Tools available

Tips & tricks

- Recommend bags with side filters to avoid inhibition.
- Draw a plate setup prior to the experiment (set lysis & PCR plates in the same order).
- Positive enrichment control handling requires separate equipment and working place.
- Work on a single strip at a time. Only move to next one after the first one is closed.
- The final lysis buffer can be stored for 2 weeks at 4°C ± 2°C.
- When transferring from lysis plate to PCR plate check the plate orientation.
- Take pipet tips in the same orientation of strips to keep track of pipetted wells.
- When labelling the strips do not write on the domed caps or optical caps.
- Save lysis sample in case there is an issue with PCR.

1

Enrichment procedure

- Enrich 25 g up to 375 g of sample or one item according to the respective protocol specified in the user guide and homogenise well.
Different samples might require different protocols – media, temperature and time. For more information, please see our certified matrix overview for BACGene kits.
- Transfer an aliquot from the enrichment culture from the filter bag into a sterile container using a bag or a serological pipette (under the laminar flow hood)



2

Lysis procedure – prepare lysis buffer

- Salmonella**
- Add 500 µL proteinase K (blue cap) to 4 mL of Lysis Buffer I (transparent)




- Listeria**
- Add 500 µL proteinase K (blue cap) and 500 µL lysozyme (orange cap) to 2.5mL of Lysis Buffer L (transparent)



3

Dispense lysis buffer - 1 strip at a time
(Leave well A1 + B1 free for negative and positive controls)

- Salmonella**
- Pipette 90 µL of Lysis Buffer I with proteinase K to each well of the lysis plate
 - Transfer 10 µL of enrichment from sterile container into the appropriate wells
 - Seal the strips with domed caps using the capping / uncapping tool
- Listeria**
- Pipette 70 µL of Lysis Buffer L with proteinase K and lysozyme to each well of the lysis plate
 - Transfer 30 µL of enrichment from sterile container into the appropriate wells
 - Seal the strips with domed caps using the capping / uncapping tool



BACGene *Listeria* monocytogenes, *Listeria* spp., *Listeria* Multiplex, *Salmonella* spp.

4

Thermal / Enzymatic lysis



- Place lysis strips/plate into 37°C (± 2°C) heating block for 20 min.



- Place lysis strips/plate into 95°C (± 2°C) heating block for 10 min.



- Transfer lysis strips/plate into pre-chilled 96 well cooling block (5°C ± 3) for 5 min.



- Spin down lysis strips/plate using a 96-well plate centrifuge or a centrifuge with adapter for 0.2 mL in 8-strip format for 30 s at 400 - 2000 x g. If condensate remains in the lids, please repeat this step.

5

PCR



- Spin down PCR strips.
- Open one strip at a time from the lysis plate and the PCR plate.



- Transfer 5 µL of sample and controls into the PCR strips.
- Spin down PCR strips.



- Transfer PCR strip to real-time instrument and start the run.
- Review results according to user guide instructions.

Before using the kit please read the instruction for use carefully for detailed instructions.

For any question, please contact your personal application specialist.

BACGene *Listeria* monocytogenes, *Listeria* spp., *Listeria* Multiplex, *Salmonella* spp.



Technologies

BACGene Protocol Overview
 certified by NF VALIDATION (AFNOR)
 and AOAC® PTMSM



BACGene Salmonella spp.
 (Cat. no. 5123221801)

BACGene Listeria spp.
 (Cat. no. 5123222101)

BACGene Listeria monocytogenes
 (Cat. no. 5123222001)

BACGene Listeria Multiplex
 (Cat. no. 5123221901)




		BACGene Salmonella spp. - AFNOR Certified Protocols						
Sample	Matrix	Milk and dairy (25 g)	All foods (25 g)	Feed (25 g)	Pet food (up to 375 g)	Milk powders & infant formula (up to 375 g)	Milk powders & infant formula + probiotics (up to 375 g)	Environmentals
Enrichment	Ratio	1:10 BPW			1:10 BPW pre-warmed		1:10 BPW ² pre-warmed	Swab (10 mL)/sponge (100 mL)/wipe/cloth (225 mL)/dust, process water (1:10) - BPW
	Temp. °C	41,5 ± 1°C			37 ± 1°C			
Time	Time	18 - 24 h	16 - 24 h	18 - 24 h				
	Prep	[Proteinase K] into Lysis Buffer I						
Lysis	Ratio	10 µL Enrichment + 90 µL Final Lysis Buffer I						
	Temp. °C	20 Minutes @37 ± 2°C + 10 Minutes @95 ± 2°C						
Plating	Ratio	5 µL Lysate on Ready-To-Use plates						
	Time	100 Minutes						



		BACGene Listeria spp., monocytogenes & Multiplex - AFNOR Certified Protocols		
Sample	Matrix	All foods (25 g)	Dust & process water (25 g/mL)	Environmentals
Enrichment	Ratio	1:10 Actero™ pre-warmed		Swab (10 mL)/Sponge (100 mL)/wipe/cloth (225 mL) Actero™ pre-warmed
	Temp. °C	37 ± 1°C		
Time	Time	18 - 24 h		
	Prep	[Proteinase K + Lysozyme] into Lysis Buffer L		
Lysis	Ratio	30 µL Enrichment + 70 µL Final Lysis Buffer L		
	Temp. °C	20 Minutes @37 ± 2°C + 10 Minutes @95 ± 2°C		
Plating	Ratio	5 µL Lysate on Ready-To-Use plates		
	Time	100 Minutes		

PCR **BAC** *Gene*
Technical

BAC Gene *Salmonella* spp.:



- The amplified target fragment is:
 - ✓ Detected with **FAM**™ fluorescence-labelled hybridization probe
 - ✓ Quenched by non-fluorescent Tide Quencher™ 2 (TQ2)
- An internal positive control (IPC) is included in the MasterMix and amplified in parallel:
 - ✓ Detected with **Cy5**™ fluorescence-labelled hybridisation probe
 - ✓ Quenched by non-fluorescent Tide Quencher™ 3 (TQ3)
 - ✓ IPC detection indicates the proper functioning of the PCR



<https://www.biosearchtech.com/gpcr-multiplex-spectral-overlay-tool>

- For every PCR, it is necessary to prepare a positive (C+) and a negative (C-) control reaction.
- Preparation of a positive (E+) and a negative (E-) enrichment controls are recommended.
- Different BACGene detection kits can be analyzed in one run

Example of simultaneous RT-PCR run for *BACGene Salmonella spp* and *BACGene Listeria Monocytogenes*:

	1	2	3	4	5	6	7	8	9	10	11	12
A	C+	7	15							C+	7	15
B	C-	8	16		BACGene <i>Salmonella spp</i>					C-	8	16
C	1	9	17							1	9	17
D	2	10	18							2	10	18
E	3	11	19							3	11	19
F	4	12	20		BACGene <i>Listeria Monocytogenes</i>					4	12	20
G	5	13	E-							5	13	E-
H	6	14	E+							6	14	E+

C+ = positive control; C- = negative control; E- = negative enrichment control;
E+ = positive enrichment control

Kit components: BACGene Salmonella spp.

Cat. no. 5123221801 of **BACGene Salmonella spp.** (96 reactions)


For Lysis:

1x **Lysis plate** for sample preparation, empty, rippable (high profile)

1x **Domed caps**, for use with Lysis plate, set of 12 strips

2x **Lysis buffer I**, 4 mL, store at $-20^{\circ}\text{C} \pm 2^{\circ}\text{C}$

2x **Proteinase K**, vials with blue cap, each with 500 μL , store at $-20^{\circ}\text{C} \pm 2^{\circ}\text{C}$




For PCR:

1x **BACGene Salmonella spp. PCR plate**, with pre-dispensed MasterMix and PCR support grid. Store light protected at $-20^{\circ}\text{C} \pm 2^{\circ}\text{C}$

1x **Optical caps**, for use with PCR plate (1 bag with 120 strips)

2x **Salmonella positive control plasmid DNA**, vial with yellow cap, 50 μL , store at $-20^{\circ}\text{C} \pm 2^{\circ}\text{C}$. Do not freeze/thaw more than 6 times

Overview: BACGene Salmonella spp.



- **Certified assays** by independent organisations (**AFNOR, AOAC**)
- **No spectral overlapping**: Spectra from Salmonella target sequence (FAM) and the internal control IPC (CY5) do not overlap
- **No internal passive reference** necessary (ROX)
- **Enriched samples** can be stored up to **72 h at 4 °C**
- **Simple lysis**: Thermal enzymatic incubation
- **Final lysis buffer** can be stored for **2 weeks at 4 °C or 2 months at -20 °C**
- **Easy handling**: Very few pipetting steps
- **Prefilled RT-PCR plate** with **ready-to-use Mastermix where IPC is included**
- **Rippable strips**; individual wells can be used
- **Simultaneous RT-PCR run**: Different pathogens can be analyzed in one run
- **Ready to use RT-PCR run templates**
- **Simple evaluation**: possible on different computers
- **Colored evaluation result display**



Q&A