ARGENTA

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Inherit DIRUI Classical Principle

- Real images - Trustworthy

No gather
No overlap

Sheath flow technology
+
High Speed Photography Technology
+
Artificial Intelligence Identification Technology (AII)

Artificial Intelligent Identification Technique

Inherit DIRUI Classical Principle

The prompt diagnosis information can be shown more directly by the report with pictures of pathological items.

- RBC:** Identify RBC subclasses automatically to indicate the source of hematuria effectively.

MIRBC NRBC ARBC SRBC OSRBC

- Cast:** In addition to identify hyaline cast, 3 pathological casts including granular cast, waxy cast and broad cast can be identified automatically to indicate the degree of kidney damage.


HYAL GRAN BROAD WAXY OCAS

- Crystal:** To indicate pyelolithiasis, renal calculi, ureterolithiasis, urolithiasis, vesical calculi, urethral calculi and liver disease.

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





All page Current page 1

Jump Previous category Next category Previous 1 totally Next Undo Back

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
Products Type	Test Item
FUS-10 II	urobilinogen, bilirubin, ketone, blood, protein, nitrite, leukocyte, glucose, specific gravity, pH
FUS-11 II	urobilinogen, bilirubin, ketone, blood, protein, nitrite, leukocyte, glucose, specific gravity, pH, ascorbic acid
FUS-11MA II	urobilinogen, bilirubin, ketone, blood, protein, pH, nitrite, leukocyte, glucose, specific gravity, microalbumin
FUS-12MA II	urobilinogen, bilirubin, ketone, blood, protein, nitrite, leukocyte, glucose, specific gravity, pH, microalbumin, ascorbic acid
FUS-13Cr II	urobilinogen, bilirubin, ketone, blood, protein, pH, specific gravity, nitrite, glucose, leukocyte, ascorbic acid, microalbumin, creatinine
FUS-14Ca II	urobilinogen, bilirubin, ketone, blood, protein, pH, specific gravity, nitrite, glucose, leukocyte, ascorbic acid, microalbumin, creatinine, Ca

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**FUS-II
Urinalysis Strips**

FUS-14Ca II
For In Vitro Diagnostic Use

Test item:
urobilinogen, bilirubin, ketone, creatinine, blood, protein, microalbumin, nitrite, leukocytes, glucose, specific gravity, pH, ascorbic acid, Ca.

Important:
1. Keep away from light and moisture.
2. Promptly replace cap after taking out strips.
3. Do not remove desiccants.

REF 3008561

LOT

 DIRUI Industrial Co., Ltd.
95 Yunhe Street, New & High Tech, Development Zone
Changchun, Jilin 130012 P.R.China
Tel: +86(431) 85100409
Fax: +86(431) 85172581
E-mail: dirui@dirui.com.cn
Http://www.dirui.com.cn

100 Strips

UROBILINOGEN 60 Sec.	3.4	normal	17	34	68	135	umol/L	
BILIRUBIN 60 Sec.	neg.		17	51	103		umol/L	
KETONE 60 Sec.	neg.	0.5	1.5	3.9	7.8	16	mmol/L	
CREATININE 60 Sec.	0.9	4.4	8.8	17.7	26.5		mmol/L	
BLOOD 60 Sec.	neg.	non hemolyzed ca.10	hemolyzed ca.10	ca.25	ca.80	ca.200	cells/ μ L	
PROTEIN 60 Sec.	10	Trace	0.3	1.0	3.0	\geq 20.0	g/L	
MICROALBUMIN 60 Sec.	neg.			0.125	0.25		mg/dL	
NITRITE 60 Sec.	neg.		ca.15	ca.70	ca.125	ca.500	cells/ μ L	
LEUKOCYTES 60 Sec.	neg.		2.6	5.6	14	26	56	mmol/L
GLUCOSE 60 Sec.	1.000	1.005	1.010	1.015	1.020	1.025	1.030	
SPECIFIC GRAVITY 60 Sec.	5.0	6.0	6.5	7.0	7.5	8.0	8.5	
pH 60 Sec.	0		0.6	1.4	2.8	5.7	mmol/L	
ASCORBIC ACID 60 Sec.	\leq 1.0		2.5	5.0	7.5	\geq 10	mmol/L	
CA 60 Sec.								
CALIBRATION AREA								


Rev.05/2019

Please read package insert carefully before use

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