

Features

- Ease of use
- Expansive and unique menu
- Accuracy comparable to reference laboratories
- Intelligent Quality Control (iQC)
- Little to no maintenance
- Detects sample interference - hemolysis, lipemia, and icterus
- Small footprint and portable (only 11.6 lbs.)
- Cost effective
- Speed: 12 minutes to results
- Stores 5,000 patient records
- Connectivity options

Benefits

- Reduce training time and increase employee compliance to use
- Save time, counter space, and costs with ability to test chemistry and electrolyte analytes on a single panel and have access to specialty tests including T4, Bile Acids, TCO₂, and Canine Heartworm on one instrument
- Diagnose and treat patients confidently
- No more worrying about running controls and calibrations as iQC monitors the analyzer, rotor, and sample before and during every analysis to ensure reliable results
- Save time and costs plus increase employee compliance to use
- Diagnose and treat patients confidently as analyte results are suppressed if affected by > 10% interference and receive useful information about sample integrity and handling
- Use this analyzer anywhere and anytime in various environments
- Save money with low-price reagents, no control and calibration costs, and no service contracts
- Obtain results, diagnose, discuss with clients and begin treatment within one visit
- Provides peace of mind that results are secure and allows recall of patient history
- Connect with Practice Management Software and external printers using the four USB ports and LAN port

VetScan VS2

Easy as 1-2-3

For more information on the VetScan VS2, contact your Abaxis Representative at 800.822.2947, email us at vetscan@abaxis.com or click www.abaxis.com/veterinary.

Abaxis is your total laboratory partner, combining the best in on-site diagnostics with expert reference laboratory services. Abaxis offers a single, comprehensive and trusted partner for all of your testing needs. The VetScan name has become synonymous with convenience, ease of use and cost-effectiveness, delivering on our promise of speed, accuracy and reliability. It's not just better technology. It's a better way.

For questions:

[800.822.2947](tel:800.822.2947) www.abaxis.com/veterinary vetscan@abaxis.com



Abaxis, Inc.
Worldwide Headquarters


3240 Whipple Road
Union City, CA 94587
Tel 800.822.2947
Fax 510.441.6150
Email vetscan@abaxis.com

www.abaxis.com

ABAXIS Europe GmbH

Bunsenstr. 9-11
64347 Griesheim Germany
Tel +49 6155 780 210
Fax +49 6155 780 2111
Email abaxis@abaxis.de

VetScan
VS2 Chemistry Analyzer



BETTER + BETTER

BETTER + BETTER.




On-Site Diagnostic Testing, with 10 Panels to Choose From

The VetScan VS2 is a state-of-the-art chemistry, electrolyte, immunoassay and blood gas analyzer that delivers uncompromising accuracy from just two drops of whole blood, serum or plasma. The VetScan VS2 is amazingly simple and intuitive, with precision reference laboratory quality results in 12 minutes—making it ideal for veterinary clinics, mobile practitioners, research laboratories, and pharmaceutical and biotech companies.

Comprehensive Test Portfolio – 10 Different Profiles to Choose From

Analytes	ALB	ALP	ALT	AMY	AST	BA	BUN	Ca	CHOL	CHW	CK	CL-	CRE	GGT	GLOB*	GLU	K+	Mg	Na+	PHOS	T4	TBIL	tCO ₂	TP	UA
Canine Wellness Profile** #500-0044	•	•	•				•	•		•			•		•	•				•		•		•	
Comprehensive Diagnostic Profile #500-0038	•	•	•	•			•	•					•		•	•	•			•	•		•	•	
Prep Profile II #500-0026		•	•				•						•			•								•	
Critical Care Plus #500-0042			•				•					•	•			•	•			•			•		
T4/Cholesterol Profile #500-0037									•													•			
Kidney Profile Plus** #500-0045	•							•	•			•	•			•	•			•	•		•		
Mammalian Liver Profile #500-0040	•	•	•			•	•		•					•									•		
Equine Profile Plus #500-0043	•				•		•	•			•		•	•	•	•	•			•			•	•	•
Avian/Reptilian Profile Plus #500-0041	•				•	•		•			•					•	•			•	•		•	•	
Large Animal Profile #500-0023	•	•			•		•	•			•			•						•			•	•	

*Calculated value **For VetScan VS2 Only



Intelligent Quality Control (iQC)

With the VetScan VS2 Chemistry Analyzer, quality control is built-in—thanks to a sophisticated iQC system that combines with proprietary algorithms to assure quality dependability with every result. No more worrying about running controls and calibrations.



Simple Operation

Practices everywhere have discovered how the VetScan Chemistry Analyzer has improved practice efficiencies and saved them time and money. Simplicity, speed, through-put with results on-site, in minutes.



1 Insert sample.



2 Insert rotor.



3 Read results.

Requires only two drops of whole blood.

