

# AUTOMATED TUBES AND BOTTLES FILLER

Versatility, accuracy and speed

# POLYWEL 3D<sup>®</sup>



1000 tubes per hour / per run\*



8 racks of 120 tubes or 4 racks of 250 tubes\*

## 3D driven

Flexible and adaptable to your different containers


## Distribution nozzles :

Ø 1mm - Ø 3mm - Ø 6mm

ACCURACY → SPEED

## Autoclavable racks

To be loaded directly from its original package



End of the cycle light indicator

## Up to 20 programs, quick to start

- Accuracy from 5ml : ± 0.1 ml
- Calibration at the beginning of the run
- Calibration during the run

## Easy cleaning

stainless steel surface



\* Ø 18mm tubes (9 ml)



**Delivered with the POLYWEL 3D :** User Manual, power chord, Polypump 3D (POLW2000), Set of tubing 3.2 mm with dispensing nozzle (POLW2008), Dispensing nozzle holder (POLW2013), rack 250.



## TECHNICAL DATA SHEET

Distribution volume	From 1 to 9999,9ml
Tubes per hour (Ø 18mm)	1000
Tubes per rack (Ø 18mm)	Large rack : 250 / Small rack : 120
Programs	20
Tubing diameters	1.6 mm / 3.2 mm / 6.4 mm
Adjustable speed	✓
Calibration	Start and during the cycle
Maximum size of containers	Height : between 115 and 250mm Width : up to 190mm Depth : up to 390mm
Dimensions	1130 x 600 x 610 mm
Power supply	120 – 240V (±10%) 50 – 60 Hz / 160W
Warranty	1 year
Country of origin	European Union - France



# POLYWEL 3D®

Distribution of media, diluent or other liquids in tubes, bottles, Petri dishes, square dishes ...



## PRODUCT CODE AND DESCRIPTION

<b>DSHB3001</b>	Tryptic Soy Agar (TSA) (500 gr)	<b>DSHB3021</b>	D-Cycloserine Selective Supplement (100mg) (10 vials)
<b>DSHB3002</b>	Tryptic Soy Broth (TSB) (500 gr)	<b>DSHB3022</b>	Desoxycholate Agar (500 g)
<b>DSHB3003</b>	Sabouraud Dextrose Agar (SDA) (500 gr)	<b>DSHB3023</b>	Kligler Iron Agar (500 g)
<b>DSHB3004</b>	Buffered Peptone Water pH = 7 (500 gr)	<b>DSHB3024</b>	L T100 Broth (500 g)
<b>DSHB3005</b>	Pharmacopoeia Diluent (500 gr)	<b>DSHB3025</b>	Müller-Kauffmann Tétrathionate Novobiocine Broth (MKTTn) (500 g)
<b>DSHB3007</b>	Mac Conkey Agar (500 gr)	<b>DSHB3026</b>	Mossel Broth (500 g)
<b>DSHB3008</b>	Mac Conkey Broth (500 gr)	<b>DSHB3027</b>	Man-Rogosa-Sharpe Agar (MRS) (500 g)
<b>DSHB3009</b>	R2A Agar (500 gr)	<b>DSHB3028</b>	MYP Bacillus cereus Agar (500 g)
<b>DSHB3010</b>	Cetrimide Agar (Pseudomonas Selective Agar) (500 gr)	<b>DSHB3029</b>	Bromocresol Purpule Glucose Agar (500 g)
<b>DSHB3011</b>	Xylose Lysine Desoxycholate Agar (XLD) (500 gr)	<b>DSHB3030</b>	OGYE Agar (500 g)
<b>DSHB3012</b>	Mannitol Salt Agar (500 gr)	<b>DSHB3031</b>	Listeria Oxford Agar (to be ordered with DSHB3051) (500 g)
<b>DSHB3013</b>	Nutrient Agar (500 gr)	<b>DSHB3032</b>	Listeria Palcam Agar (to be ordered with DSHB3052) (500 g)
<b>DSHB3014</b>	Rappaport-Vassiliadis Broth (RVS) (5 KG)	<b>DSHB3033</b>	Bolton Broth Campylobacter (to be ordered with DSHB3053) (500 g)
<b>DSHB3015</b>	Brilliant Green Bile Broth (BGBLB) (500 g)	<b>DSHB3034</b>	Tryptic Glucose Yeast Agar (PCA) (500 g)
<b>DSHB3016</b>	Brain Heart Infusion Agar (500 g)	<b>DSHB3035</b>	Potato Dextrose Agar (PDA) (500 g)
<b>DSHB3017</b>	Columbia Agar (500 g)	<b>DSHB3036</b>	Sabouraud Broth (500 g)
<b>DSHB3018</b>	Fraser Broth (500 g) (to be ordered with DSHB3049 or DSHB3050)	<b>DSHB3037</b>	Sabouraud Dextrose Agar (SDA) + chloramphenicol (500 g)
<b>DSHB3019</b>	Lyophilised Rabbit Plasma Fibrinogen (RPF Supplement) (4 x 100 ml)	<b>DSHB3038</b>	Tryptone Bile X-Glucuronide Agar (TBX) (500 g)
<b>DSHB3020</b>	Pseudomonas CFC Supplement (10 x 500 ml)	<b>DSHB3039</b>	TCBS Kobayashi Agar (500 g)



# MEDIAFEED®

## Dehydrated Media

### PRODUCT CODE AND DESCRIPTION

<b>DSHB3040</b>	Thioglycolate Broth (500 g)
<b>DSHB3041</b>	Peptone Salt Broth (500 g)
<b>DSHB3042</b>	Tryptone-Sulfite Cycloserine Agar (TSC) (to be ordered with DSHB3021) (500 g)
<b>DSHB3043</b>	Violet Red Bile Glucose Agar (VRBG) (500 g)
<b>DSHB3044</b>	Violet Red Bile Lactose Agar (VRBL) (500 g)
<b>DSHB3045</b>	Chloramphenicol Glucose Yeast Extract Agar (YGC) (500 g)
<b>DSHB3046</b>	Bile Esculine Azide Agar (500 g)
<b>DSHB3047</b>	Endo Agar (500 g)
<b>DSHB3048</b>	Buffered Peptone Water pH7.2 (500 g)
<b>DSHB3049</b>	Fraser Selective Supplement (10 x 500 ml)
<b>DSHB3050</b>	Fraser Half Selective Supplement (10 x 225 ml)
<b>DSHB3051</b>	Oxford Antimicrobial Supplément (10 x 500 ml)
<b>DSHB3052</b>	Palcam Antimicrobial Supplément (10 x 500 ml)
<b>DSHB3053</b>	Bolton Broth Selective Supplement (10 x 500 ml)
<b>DSHB3054</b>	Mueller Hinton (500 g)
<b>DSHB3055</b>	Rappaport Vassiliadis Soy Broth (RVS) (500 g)
<b>DSHB3056</b>	Columbia Agar Base (500 g)



**FEED YOUR MEDIAWEL®!**

