

Control to fit your amplification



Thermal Cycler Fleet Control Software

Manage multiple users, methods, and instruments with a single intuitive interface


Applied Biosystems™ Thermal Cycler Fleet Control Software is as intuitive and easy to use as our thermal cyclers. With this tool, you'll be able to view and control all Applied Biosystems™ thermal cyclers.

Secure—Thermal Cycler Fleet Control Software comes preinstalled on a server that sits securely behind your firewall

Flexible—create custom permissions and rules for users, methods, and instruments

Powerful—quickly view your thermal cycler inventory and availability right from the dashboard

Convenient—control your thermal cyclers, reports, methods, and firmware upgrades from your office, another lab computer, or any other networked device




Enjoy peace of mind knowing that Thermal Cycler Fleet Control Software operates completely inside of your organization's firewall.



Superior service and support

Help reduce downtime with rapid exchange support



Thermal Cyclers Fleet Control Software arrives preinstalled on a server and can be operated from the server or from any computer or device on the same network using the Google™ Chrome™ internet browser.

Specifications

Compatible Applied Biosystems™ instruments

ProFlex™ PCR System; VeritiPro™ Thermal Cyclers; SimpliAmp™ Thermal Cyclers; Automated Thermal Cyclers; and MiniAmp™ Thermal Cyclers

Key capabilities

- View all instruments together on a dashboard
- Automatically store and easily access run reports
- Manage multiple user accounts with password protection
- Track system changes with an audit log
- Identify methods or instruments via barcode data entry

Hardware

Dell™ OptiPlex™ 3050 All-in-One Computer server with touchscreen, keyboard, and mouse

Hardware memory

128 GB solid-state hard drive

Hardware display interface

19.5-inch wide-viewing angle, HD, 1,600 x 900 resolution

Hardware power

External power supply unit; 130 W adapter

Network connections

Ethernet LAN 10/100/1,000 or Wi-Fi

Compatible browser

Chrome

Standard warranty

3 years

Ordering information

Product

Cat. No.

Thermal Cyclers Fleet Control Software

A40070

ProFlex PCR System, 96-well*

4484075

SimpliAmp Thermal Cyclers*

A24811

VeritiPro Thermal Cyclers, 96-well*

A48141

Automated Thermal Cyclers, 96-well*

A31486

MiniAmp Thermal Cyclers*

A37834

* Thermal Cyclers Fleet Control Software is compatible with all versions of Applied Biosystems thermal cyclers.

Find out more at thermofisher.com/fleetcontrol

ThermoFisher
SCIENTIFIC

For Research Use Only. Not for use in diagnostic procedures. © 2021 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. Dell and OptiPlex are trademarks of Dell Inc. Google and Chrome are trademarks of Google Inc. **COL016925 0921**