

## Thermo Scientific Digital Tube Revolver

Increase your lab's efficiency with our variable speed tube revolver. Its small footprint, flexible design, and simple display help increase flexibility and effectiveness.



- Interchangeable rotisseries accommodate different tube sizes for maximum flexibility. Each paddle can be adjusted from 0-90 degrees.
- Achieve more effective mixing with “agitation” mode. Agitation will occur at every 40 degrees in the rotation cycle.
- Small footprint
- Plug-and-play simplicity
- Maintenance-free motor, quiet operation

### Common Applications

- Molecular Biology
- Clinical
- Biochemistry



Cat. No.	88881001	88881002
<b>Plug Configurations</b>		
<b>Rotation Speed</b>	10-40 rpm	
<b>Capacity</b>	Combination of 4 paddles of: 7 x 10/15ml tubes 7 x 5/7ml tubes 15 x 1.5ml/2ml tubes 21 x 0.5/0.8ml tubes One paddle of 6 x 50ml tubes	
<b>Standard Rotisserie Paddles included</b>	<b>88881005</b> – 10/15ml paddle, quantity 2 <b>88881006</b> – 5ml/7ml paddle, quantity 2 <b>88881007</b> – 1.5ml/2ml paddle, quantity 2 <b>88881008</b> – 0.5ml/0.8ml paddle, quantity 2 <b>88881009</b> – 50ml paddle, quantity 1  Paddles can also be purchased as accessories in quantities listed above	
<b>Operating Temperature</b>	4°C – 60°C	
<b>Dimensions W x L x H</b>	225 x 154 x 195mm / 8.9 x 6 x 7.7in	
<b>Max Power</b>	6W	
<b>Weight</b>	1.2kg / 2.65lbs	
<b>Power</b>	AC100V-240V / 50Hz – 60Hz	
<b>Certification</b>	cCSAus, CE	
<b>Warranty</b>	2 years	

Ordering Information and Specifications	Cat. No.
Paddles 10ml/15ml, 2 pack	<b>88881005</b>
Paddles 5ml/7ml, 2 pack	<b>88881006</b>
Paddles 1.5ml/2ml, 2 pack	<b>88881007</b>
Paddles 0.5ml/0.8ml, 2 pack	<b>88881008</b>
Paddle 50ml	<b>88881009</b>

