



ARGENTA
WET

LOCMEDT®

Automatyczny analizator biochemiczny

Noahcali-100



W pełni automatyczny - wstępnie zainstalowany rozcieńczalnik, kalibracja na podstawie kodu QR



Analiza do 26 wskaźników biochemicznych jednocześnie



Wymagana mała objętość próbki - jedynie 90-120 μ l



Dokładne wyniki już w 12 minut



Automatyczny analizator biochemiczny

Noahcali-100 jest przenośnym analizatorem biochemicznym umożliwiającym ilościowe oznaczanie wskaźników biochemicznych z krwi pełnej, osocza i surowicy.

Jednorazowe dyski z zliofilizowanymi odczynnikami oparte są na technologii mikroprzepływowej i umożliwiają uzyskanie wyniku z 3 kropli krwi.

Trzy proste kroki analizy:



Dodaj 100 µl próbki



Włóż dysk



Odczytaj wynik

Produkty	L. szt./op.	Nr kat.
Health Check Panel 16 (TP, ALB, CK, A/G*, TBIL, ALT, BUN, CRE, GLU, TC, Ca, PHOS, GLOB*, AMY, ALP, BUN/CRE*)	10	35200016
Liver Function Panel 13 (ALB, AST, ALT, ALP, GGT, TBA, TBIL, TC, TP, UA, GLOB*, A/G*, AST/ALT*)	10	35200020
Electrolytes Panel 17 (BUN, Ca, P, Mg, Na+, K+, Cl-, GLU, BUN/CRE*, Na/K*, CREA, tCO2, LAC, PH, AG*, Ca*, PHOS*, Osm*)	10	35200076
Kidney Function Panel 9 (ALB, PHOS, BUN, CRE, UA, Ca, K+, tCO2, BUN/CRE*)	10	35200021
Liver and Kidney Function Panel 14 (ALB, AST, ALT, ALP, CRE, GLU, GGT, TBIL, TP, BUN, GLOB*, A/G*, BUN/CRE*, AST/ALT*)	10	35200019
Peanesthetic Panel 14 (AST, ALT, ALP, CK, CRE, GLU, LDH, TP, BUN, BUN/CRE*, AST/ALT*, ALB, GLOB*, A/G*)	10	35200017
Comprehensive Diagnostic Panel 26 (TP, ALB, GLOB*, A/G*, TBIL, ALT, AST, BUN, CRE, TBA, GLU, TG, TC, CA, GGT, LDH, PHOS, Ca*, P, CK, AMY, ALP, BUN/CRE*, AST/ALT*, LPS, UA, C-Ca*)	10	35200035
Critical Care Panel 12 (ALT, AMY, CHE, BUN, CRE, GLU, NA+, K+, tCO2, BUN/CRE*, Na+/K+*, CL-)	10	35200018
Large Animal Panel 21 (ALB, AST, ALP, CK, CRE, PHOS, LDH, GGT, Ca, TBIL, GLU, TG, TP, BUN, Na+, K+, Mg2+, GLOB*, BUN/CRE*, A/G*, Na+/K+*)	10	35200027
Health Check Plus Panel 20 (ALB, ALT, AMY, ALP, CRE, PHOS, TC, GGT, Ca, TBIL, GLU, TP, BUN, Na+, K+, tCO2, GLOB*, BUN/CRE*, A/G*, Na+/K+*)	10	35200028
Coagulation Panel 4 (APTT, PT, TT, FIB)	10	35200031
Canine Inflammation Panel 10 (CRP, LPS, AMY, TP, ALB, A/G*, GLO*, BUN, CRE, BUN/CRE*)	10	35200039
Feline Inflammation Panel 13 (SAA, LPS, AMY, ALB, TP, A/G*, GLO*, BUN, CRE, BUN/CRE*, TBA, GGT, ALP)	10	35200060



ARGENTA

Argenta Sp. z o.o., ul. Człuchowska 6, 60-434 Poznań
t. +48 61 847 46 37 e. info@argenta.com.pl w. www.argenta.com.pl