



ARGENTA

Chemia kliniczna - parametry

BS-600M, BS-240 PRO

mindray





Panel wątrobowy

Test	Pełna nazwa	Metoda	Test	Pełna nazwa	Metoda
AST/GOT	Aspartate Aminotransferase	IFCC Method	Bil-T(VOX)	Bilirubin Total	VOX Method
ALT/GPT	Alanine Aminotransferase	IFCC Method	TP	Total protein	Biuret Method
ALP	Alkaline Phosphatase	IFCC Modified Method	TP II	Total protein II	Biuret Method
GGT	Gamma- Glutamyltransferase	Szasz Method/IFCC stand	CHE	Cholinesterase	DGKC Method
ALB	Albumin	Bromcresol Green Method	PA	Prealbumin	Turbidimetry Method
Bil-D(DSA)	Bilirubin Direct	DSA Method	TBA	Total Bile Acids	Enzymatic cycling assay
Bil-D(VOX)	Bilirubin Direct	VOX Method	AFU	α -L-Fucosidase	CNPF Method
Bil-T(DSA)	Bil-T(DSA)	DSA Method			

Panel nerkowy

Test	Pełna nazwa	Metoda
UREA	Urea	Urease-GLDH, UV Method
UA	Uric acid	Uricase-Peroxidase Method
CREA-S	Creatinine	Sarcosine Oxidase Method
MALB	Microalbumin	Immunoturbidimetric Assay Method
TPUC	Total Protein in Urine/CSF	Pyrogallol Red- Molybdate Method
Cys-C II	Cystatin C II	Latex Immunoturbidimetric Method
RBP	Retinol Binding Protein	Particle-enhanced Immunoturbidimetric Assay Method
B2-mG II	B2- Microglobulin	Latex Immunoturbidimetric etric Method

Panel zapalny

Test	Pełna nazwa	Metoda
CRP	C-Reactive Protein	Turbidimetry Method
ASO II	Antistreptolysin "O"	Latex Immunoturbidimetric etric Method
RF II	Rheumatoid Factor	Immunoturbidimetric etric Method
IgA	Immunoglobulin G	Turbidimetry Method
IgG	Immunoglobulin G	Turbidimetry Method
IgM	Immunoglobulin M	Turbidimetry Method
IgE	Immunoglobulin E	Particle-enhanced Immunoturbidimetric etric Assay Method
C3	Complement C3	Turbidimetry Method
C4	Complement C4	Turbidimetry Method

Panel anemii

Test	Pełna nazwa	Metoda
FER	Ferritin	Particle-enhanced Immunoturbidimetric etric Assay Method
TRF	Transferrin	Immunoturbidimetric etric Assay Method
UIBC	Unsaturated Iron Binding Capacity	Colorimetric Method
G6PD	Glucose-6- Phosphate Dehydrogenase	UV Enzymatic Method



Panel lipidowy

Test	Pełna nazwa	Metoda
TC	Total Cholesterol	CHOD-PO Method
TG	Triglycerides	GPO-POD Method
HDL-C	HDL- Cholesterol	Direct Method
LDL-C	LDL- Cholesterol	Direct Method
Apo-A1	Apolipoprotein A1	Turbidimetry Method
Apo-B	Apolipoprotein B	Turbidimetry Method
Lp(a) II	Lipoprotein (a)	Latex Immunoturbidimetric etric Method

Panel sercowy

Test	Pełna nazwa	Metoda
CK	Creatine Kinase	IFCC Method
CK-MB	Creatine Kinase-MB	IFCC Method
LDH	Lactate Dehydrogenase	IFCC Method
α -HBDH	α -Hydroxybutyrate Dehydrogenase	DGKC Method
HS-CRP	High Sensitivity C-reaction Protein	Particle-enhanced Immunoturbidimetric
MYO	Myoglobin	Particle-enhanced Immunoturbidimetric Assay Method
D-Dimer	D-Dimer	Particle-enhanced Immunoturbidimetric Assay Method
HCY	Homocysteine	Enzymatic Cycling Method

Panel cukrzycowy

Test	Pełna nazwa	Metoda
Glu-GOD	Glucose	GOD-POD Method
Glu-HK	Glucose	HK Method
FUN	Fructosamine	Colorimetric Assay
HbA1c	Hemoglobin A1c	Enzymatic Assay Method
B-HB	B- Hydroxybutyrate	Enzymatic Colorimetric Method

Panel elektrolitów

Test	Pełna nazwa	Metoda
Mg	Magnesium	Xylidyl Blue Method
P	Phosphorus	Phosphomolybdate Method
P II	Phosphorus II	Phosphomolybdate Method
Ca	Calcium	Arsenazo Method III
Fe	Iron	Colorimetric Assay
CO2	Carbon dioxide	Enzymatic method

Panel płucny

Test	Pełna nazwa	Metoda
ACE	Angiotensin Converting Enzyme	Enzymatic Colorimetric Assay Method
ADA	Adenosine Deaminase	Enzymatic Colorimetric Method

Panel trzustkowy

Test	Pełna nazwa	Metoda
α -AMY	α -Amylase	IFCC Method
LIP	Lipase	Enzymatic Colorimetric Assay Method

